

## UROSKOP D

**SP**

### Service Instructions

Standard table pads  
for metal table extensions

© Siemens AG 1997

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

Print No.: RLL5-310.076.08.01.02

Replaces: n.a.

English

Doc. Gen. Date: 02.97

## General

- Table pads, Item nos. 3750 366 (27 cm) and 3750 374 (83 cm), can be attached to the older tabletop extensions (metal construction ) using Velcro pads.

### NOTICE

Table pads can be attached to the carbon fiber table extensions Item nos. 1676 910 (83 cm) and 1676 928 (27 cm) as well as metal table extensions using this procedure.

## Material

- 4 pieces (fuzzy) Velcro pads, Item no. 3838 344 per table pad

## Service equipment

- Measuring tape.

## Procedure

- To determine the position for the adhesive:
  - Move the table in the horizontal position and attach the tabletop extension.
  - Center the pads on the table extension.
  - Mark the locations for the adhesive that correspond to the Velcro strips attached to the underside of the tablepad (refer to the drawing).
- Clean the positions where adhesive will be applied.
- Remove the protective backing from the (fuzzy) Velcro pads and attach them to the marked positions on the upper side of the table extension pressing them down firmly.

